

VITREOUS CUTTER



Cod. YD7C



- YD7C can provide a rapid and stable response of vacuum by using a venture pump as vacuum power.
- Top cut frequency can attach 2500 CPM, and at the same time, it can carry out micro-cutter by reducing the pulling of around tissues.
- Audio system: according to the cutter frequency, the intensity of sound is proportion to vacuum.
- Online pipe sensing technology can provide a stable vacuum for operator through controlling system.
- 12.1 inch LCD touch screen can make vacuum of irrigation/aspiration(I/A), air-liquid exchange pressure, out-put diathermy power or volume of cutter frequency etc. displayed in real time.

Technical features	
Mainframe size	530 x 410 x 650 mm
Mainframe weight	32kg
Footswitch control size	400 x 200 x 140 mm
Pushcart size	450 x 400 x 700 mm
Pushcart weight	7.5kg
Application voltage	(use switch to exchange) 220±22V, 50±1Hz/110±11V, 60±1Hz
Input power	≤ 250VA
Fuse	choose according to the power 250V/3.15A(51S-032H) or 125V/5A(51S- 050L)
Pneumatic pressure	0.5-0.8Mpa, Nitrogen or compressed air(no use other kinds of gases)
Modes of Vitreous cut	Scope of cut frequency is from 60-2500 CPM, contains 3 continuously adjustable forwards, which are 60-800, 850-1600, 1650-2500 CPM, with instant switching between three preset cut rates.
Vacuum I/A mode	presetting vacuum volume is 133- 798hPa(100-600mmHg), contains 50 adjustable forwards, controlled by footpedal with vacuum volume.
Air-liquid exchange mode	preset 19 adjustable forwards 13.3hPa~133hPa (10-100mmHg)
Plumb Gas Scissor Model	both proportion and single scissor are optional
Working Model of Bipolar Diathermy	$300 kHz$, 600Ω load, Maximal output power is $5W$, and the output power can be adjusted for $5\% \approx 100\%$
Audio Frequency Clews	Vacuum Clew, Diathermy Clew, Warning Clew of gas-liquid exchange vacuum.

Pulse Healthcare Technology House: 11/1, Shahid Minar Road Kallyanpur, Dhaka-1216, Bangladesh. Mobile: +8801708008061

Mobile: +8801708008061 e-mail: info@pht.com.bd web: www.pht.com.bd



ESSE3 srl, Via Garibaldi 30 14022 Castelnuovo D.B. (AT) Tel +39 011 99 27 706 Fax +39 011 99 27 506 e-mail esse3@chierinet.it web: www.esse3-medical.com

