

YZ20T4 Operating Microscope





YZ20T4



Operating Microscope



- Assistant microscope can be rotated to opponent microscope.
- Internal beam splitter can be connected with teaching tube and CCD system.
- Equipped with three different objective lenses (f200/f250/f300) and applicable in ophthalmic surgery, orthopedic surgery and hand surgery.
- Oblique illumination can also be used as slit illumination and slit width is adjustable

Technical Features	
- M 10 11	12.5 (100)
Eyepiece Magnifications	12.5×/18B
Diopter Adjustment	±7D
Pupil Distance	50~75mm
Magnifications for Main Microscope	f200: 4X 6X 10X 16X 25X
	f250: 3.2X 5X 8X 12.8X 20X
	f300:2.7X 4.2X 6.7X 11X 17X
Magnifications for Assitant Microscope	f200: 6X 10X 16X
	f250: 5X 8X 12.8X
	f300: 4.2X 6.7X 11X
Maximal Resolution	112 LP/mm
Illumination Source	12V/100W , Cold Reflection Halogen Bulb for Medical Use
Illumination Type	6°+0°coaxial illumination of cold light source 25°Oblique Illumination (It can be used as slit illumination) Slit Width is adjustable, and its angle is 360°rotation.
Coaxial Illumination	≥ 50000lx
Oblique Illumination	≥ 45000lx
Reaching Radius of Arm	1040mm
Adjustable Vertical Range	850mm~1350mm
Fine Focusing Speed and Range	≤2mm/s, ≥ 50mm
Speed for X/Y and Range	≤2mm/s, 50mm x 50mm
Input Voltage	AC220V±22V/50Hz±1Hz, AC110V±11V/60Hz±1Hz
Power	330VA
Electrical Safety Standard	Executive Standard: GB9706.1-2007, type I
Total Weight OISTRIBUTOR	215kg





House: 11/1, Shahid Minar Road Kallyanpur, Dhaka-1216, Bangladesh.

Mobile: +8801708008061 e-mail: info@pht.com.bd web: www.pht.com.bd



ESSE3 srl, Via Garibaldi 30 14022 Castelnuovo D.B. (AT) Tel +39 011 99 27 706 e-mail esse3@chierinet.it web: www.esse3-medical.com



