

ES-AUN Line

Bench Top Steam Sterilizers





ES-AUN Line



Bench Top Steam Sterilizers



Features

- · Automatic and safety door lock system
- LCD display, easy for operation
- Steam compensation techniques
- LCD display. 2 programs, easy to choose,121/134.
- Fully-Automatic feeding water, excess water and air to be drained away
- Display the failure reason on LCD screen. Easy to correct.
- Can show the reason of failure on LCD, don't need to check recording to code menu.

Safety Devices

- Machine has safety door lock system
- Over heat control: when inner temperature is over heated, power will be switch off
- Automatically and start alarm system
- Over pressure control: alarm will be announced when the machine is over pressure
- · Low water level detected: alarm will be announced when it is low water

Optional basket tray

Technical Features

	ES-AU14N	ES-AU17N	ES-AU23N	ES-AU30N	ES-AU45N
Chamber capacity	14 Liters	17 Liters	23 Liters	30 Liters	45 Liters
Chamber Dimension (ØxD)	25x28 cm	25x38 cm	25x49 cm	25x60 cm	36x62 cm
Overall Size (LxWxH)	48x48x42 cm	56x48x42 cm	69x48x42 cm	79x48x42 cm	84x66x56 cm
Heater	1200W	1200W	1500W	1600W	3500W
Net weight	31 kg	35 kg	38 kg	48 kg	72 kg
Gross weight	38 kg	40 kg	43 kg	58 kg	82 kg
Packing size	71x58x48 cm	71x58x48 cm	81x58x48 cm	91x58x48 cm	90x72x72 cm
Minimum table size	48x35 cm	48x40 cm	48x50 cm	75x48 cm	80x66 cm
Power supply	110V, 220V, ±10%			110V, 220V, ±10%	
Danie anno anadonna	Col. Class Boots de la color				

Only Class B autoclave

1 tray holder with 3 shelves

121°C/134°C Please, refer to different programs Stainless steel AISI 304



Stainless steel

instrument container

erilization temperatur

Pulse Healthcare Technology

House: 11/1, Shahid Minar Road Kallyanpur, Dhaka-1216, Bangladesh.

Mobile: +8801708008061 e-mail: info@pht.com.bd web: www.pht.com.bd



ESSE3 srl, Via Garibaldi 30 14022 Castelnuovo D.B. (AT) Tel +39 011 99 27 706 Fax +39 011 99 27 506 e-mail esse3@chierinet.it web: www.esse3.dreamgest.com



